

TATO

Thermal Ablation Treatments for Oncology



Erase **CANCER** in day hospital



Thermal Ablation Treatments for Oncology



1 day to eliminate the cancer

Thermal ablation is the permanent destruction of targeted tissues caused by a local increase of temperature.

Thermal ablation procedures are minimally invasive alternatives to surgical resection. They are performed in Day-Hospital and are associated with a significant reduction of operative mortality and hospitalization costs.

TATO is a compact, mobile and highly efficient multi-applicator thermal ablation device.

2 market-leading technologies

TATO is a microwave thermal ablation device.

Compared to the older radiofrequency technology, MWA can generate higher temperatures in less time since tissue charring does not hinder the radiation of the microwaves and it is less susceptible to the heat-sink effect of peritumoral vessels.

Significant reduction of side effects:



Comet effect




Charring



Chimney effect

optimal ablation shape



TATO

3 reasons to choose TATO

Efficiency: TATO and TATOpro applicators are designed to minimize energy losses.

The VLL™ (Very Low Loss) technology allows the applicators to operate, in absence of a cooling system, without causing skin burns.

Modularity: TATOpro applicators come in various sizes, ranging from 11G to 18G.

They are suitable both to perform small, highly precise ablations and to treat big multifocal tumors.

Safety: TATO has an automatic real-time safety control loop that stops the emission of MW power in case a state of malfunctioning of the applicator is detected.



From 0 to 7 in 10'

TATO is a multi-applicator thermal ablation system. By delivering a total MW power output ranging from 10 to 100 Watts, the device is capable of performing ablations as small as 0.5 cm, up to big 7 cm multifocal lesion, all in a single 10 minutes session.

The multi-applicator technology allows the user to set the MW generator to work with one-to-four applicators of mixed Gauges and lengths in order to plan the most accurate intervention.

The final lesion shape produced by a TATOpro applicator has an optimal sphericity, furthermore sphericity increases by employing two or more applicators on a single tumor mass.

All the applicators are simultaneously controlled by a smart Human Machine Interface, which allows the operator to pause the procedure and make adjustments of the treatment parameters in real-time.

Multi-applicator for multi-applications

The possibility of planning a highly accurate treatment of tumor masses of very different sizes and shapes opens up new perspectives in the thermal ablation field.

New applications are being studied on organs that rarely receive treatment with MW thermal ablation such as prostate, thyroid, breast and bones that will soon be able to benefit from all its advantages.




TATOpro applicators use the VLL™ (Very Low Loss) technology that intrinsically limits the overheating of the shaft, resulting in the absence of a complicate and unreliable cooling system.




As an additional safety measure, TATOpro applicators come with a small, magnetic temperature clamp: a special accessory that uses sensitive liquid crystal labels to monitor the temperature directly on the shaft.

Design



Very thin handles with ergonomic grip and color code

Precision



Definite and controllable lesions for one or multiple tumors

Modularity



Probes of any Gauge and length can be combined


Optimal Shape



0.9 sphericity index with two uncooled applicators



TATO_{pro}



The power of 4

Thanks to the FSS™ (Fast Sequential Switching) technology, built inside the automatic power distribution unit, TATO is able to handle up to four applicators with a single generator.

Interventional radiologists and surgeons will be able to operate big multifocal cancers or up to four dislocated metastasis in a single 10 minutes session.

TATOpro applicators are covered with Teflon, which makes them echogenic, self-lubricating and non-sticky.

TATOpro handles are the first one that come with a detachable and reusable cable, allowing an easier than ever placement of the applicator, which is especially useful for CT-scan positioning.

TATOpro applicators can be equipped with a luer lock connector for injecting drugs, alcohol, saline solution or microbubbles during the ablation treatment.

Finally, TATOpro is the first MW hollow probe that allows the insertion of an electromagnetic navigation sensor, to achieve a very high positioning accuracy.

Size



TATOpro applicators are available in 11, 14, 17 and 18G

Reusable Cable



The low loss MW cable can be detached and sterilized

Liquid Injections



TATOpro applicator allow injecting liquids into the diseased tissue

Multi Applicator



4 applicators controlled with one single interface

Single applicator treatment

A single TATOpro applicator can perform ablations of different sizes depending on three treatment parameters: Gauge, power delivered and time.



A 14G uncooled probe will produce a 3.5 x 5 cm ablation with 40 W of applied power in 10 minutes, suitable for medium sized liver tumor ablation.

A 17G uncooled probe will produce an average 2.3 x 4,2 cm ablation in 10 minutes with 20 W of applied power.


A 18G uncooled probe will produce an extremely precise ablation of 0.5 x 0.9 cm in a 2 minutes application with only 10 W of applied power, suitable for small nodule ablations.

Ablation shape possibilities

A spherical 6 x 6 cm ablation for the treatment of big multifocal tumors (3 probes).
The simultaneous treatment of multiple metastases (2, 3 or 4 probes).
Ablating the Parenchymal Transection Margin to reduce blood loss during a surgical liver resection (2, 3 or 4 probes).

Different configuration in shape and power

Big lesions basic shape with multiple needles



Available for liver resection





Simple




Easy interface design.
Full control of the ablation process



TATO can be completely controlled via an intuitive easy to use touchscreen Human Machine Interface. The user will be able to fully set up the device in few minutes following a simple 7 step wizard.




Time and power settings for each applicator can be checked and modified on a single step of the configuration wizard. A maximum of 100 W can be distributed among the applicators and the recommended treatment duration is 10 minutes.



The Human Machine Interface allows pausing the treatment and adjusting the MW power emitted by each applicator output depending on the feedback from imaging systems.

Simultaneous monitoring of all the applicators in use on the same display.



The power switching unit is designed to host from 1 to 4 applicators at the same time and to deliver the energy emitted by the generator to each one of them in the amount established by the operator.



TATO main unit comes equipped with USB port, a lit power button and two indicator lights: a yellow light that illuminates when the generator is running and a red one that will light in case a state of malfunctioning is detected.



The adjustable arm allows bringing the power switching unit closer to the bed of the patient, reducing encumbrance and allowing the use of shorter cables, which results in a further reduction of energy losses.



TATO main unit, integrating MW generator, power supply, system processor and touchscreen display, has a very compact design (36 cm x 31 cm x 13 cm) and weights 7.5 Kg.



Ergonomic handle allows moving the entire device easily with one hand.



A low loss MW coaxial cable connects TATO main unit to the Switch where TATOpro probes are connected, further reducing energy losses.



The cart is equipped with a cable reel and a container for manual and accessories.



5 self-locking wheels cart guarantees compliance with general requirements for basic safety and essential performance EN60601-1 3rd edition.



erased



TATO

biomedical

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

CE
0426